




**"EZ POT"
PORTABLE COMMODE**



**"EZ POT"
PORTABLE COMMODE**

(Drop It In The Pot)



The EZ Pot is beneficial for care providers and patients because it **reduces injuries, patient's falls, pressure ulcers, infections, bedside commodes height concerns,** and makes the experience more **comfortable.**

**NEW AGE INNOVATIONS INC.
P.O. BOX 28555
Birmingham, Alabama 35228**

**Website:
newageinnovations.org
Phone:
205-541-0885
Fax:
205-663-1062
E-mail:
lawrencemcnair@
newageinnovations.org**

**NEW AGE INNOVATIONS
INC.**

**TOUCHING LIVES
ONE INNOVATION AT
A TIME**

The EZ Pot can also be used by patients struggling with **joint pains, illness, balance problems and surgeries.**

The EZ Pot enables care providers the opportunity to **replace existing bedside commodes and bedpans.**

**(Includes: Toilet Frame,
Cushion, Tray, Urinal, Foam
Brushes, and Cleaning Bags)**

EZ Pot

Weight: 7LBS Capacity: 750 LBS

Size: 21" x 16.5" x 3.75"

(L x W x H)

Toilet Frame

Weight: 9LBS Capacity: 650 LBS

Size: 26" x 22" (W x D)

Disposable Tray

Size: 14" x 12" (L x W)

Disposable Cushion

Size: 12.5" x 16" x 2"

(L x W x H)

Disposable Cleaning Bags

(Pack of 20)

Disposable Brushes

(Pack of 12)

*** Includes 1 Year Warranty**

*** Patent Pending**

In today's society, there are issues when caring for immobile, geriatric and bariatric patients when relieving themselves. One issue involves lifting patients, which is time consuming and may cause injury to patients and care providers. Another issues involves the concern of patient falling when attempting to relieve themselves, resulting in serious injuries. The EZ Pot consist of toilet rails, three handles, an attachable padded seat cushion, a removable tray, foam brush, and cleaning bags. The EZ Pot improves comfort and make it easier for care providers to dispose of personal waste from patients. The EZ Pot reduces skin breakdown after urination and feces voids from patients. The EZ Pot also reduces time and cleanup efforts for care providers. The EZ Pot is ideal for use in hospitals, rehabilitation centers, nursing homes, assistant living, and in-home care because it makes the process much easier and less embarrassing for patients. The EZ Pot can be used as a bedside commode or in the bed, chair, or on other flat sturdy surfaces by detaching it.

Solidly Constructed

- Composed of lightweight **Lexan plastic**.

- **Three grip handles** provide care providers the support needed to maneuver the EZ Pot underneath patients.
- Reduce **"soiled"** patients' beds, which causes **infections**.
- Designed to keep bedpans from **breaking** underneath patients.
- Allows patients to sit in a **natural position (upright) in bed**, improving their ability to eliminate.
- **Attachable seat cushion, curved corners** and **smooth surface** ensures patients' comfort.
- Clean with **antibacterial cleaning agent** by using **foam brush**, or use **cleaning bags for easy cleanup**

Cost Effective

- The **EZ Pot** holds over **2.5 quarts, 72 ounces, 2 liters, and 4.5 pints** of liquid.
- The **EZ Pot** holds over **4.5 pounds** of solids.
- Reduce spending revenue on **bariatric bedside commodes** for patients needs.
- Can be used in a **bed, chair, sofa, car or other flat surfaces**.